



Radiation Products Design, Inc.

5218 Barthel Industrial Drive

Albertville, MN 55301

www.rpdinc.com

Phone: 800-497-2071 Fax: 763-497-2295

**RPD is an
authorized distributor**

RPD Product Information

Item Number Description

300-650-BNC-M * PTW Model 30004 / 30012 - 0.6cc C/A1 Farmer Chamber

IONIZATION CHAMBER TYPE 30004

Radiation Therapy

Diagnostic Radiology

Nuclear Medicine

Radiation Protection

"C/Al-Farmer" Chamber 0.6 cm³


Description

The 0.6 cm³ ionization chamber type 30004 is the chamber for absolute dosimetry for use with therapy dosemeters where a graphite/aluminium chamber is required. This chamber is of delicate construction since it is equipped with a graphite cap. It should be handled with care. The chamber is fully guarded up to the measuring volume. The outer dimensions are fully compatible with the Farmer chambers of other manufacturers.

Nominal useful range 60 keV – 50 MeV

Technical data

Volume:	0.6 cm ³
Response:	2–10 C/Gy
Leakage:	$< 4 \cdot 10^{-14}$ A
Polarizing voltage:	max. 500 V
Cable leakage:	10^{-14} C/(Gy · m)
Wall material:	Graphite
Wall density:	1.85 g/cm ³
Wall thickness:	0.425 mm
Area density:	79 mg/cm ²
Electrode:	Al 0.6 mm, 1 mm Ø, 21.2 mm long
Range of temperature:	+10 °C ... +40 °C
Range of rel. humidity:	10 % ... 80 % (< 20 g/m ³)
Ion collection time:	300 V: 0.18 ms 400 V: 0.14 ms 500 V: 0.11 ms



- Guard ring up to measuring volume
- Suitable for use in solid state phantoms, watertight sleeves available
- Guard ring at potential of the collecting electrode
- Touchable parts free of high voltage
- High voltage to be connected only with active current limiting device ($I_{max} < 0.5$ mA)
- Open measuring volume, without check device air density correction is necessary
- Lengthening cable up to 100 m available
- Connector: TNC or BNC

Saturation behaviour

	Polarizing voltage	99.0 % saturation	99.5 % saturation
Max. dose rate at continuous irradiation	300 V	5.7 Gy/s	2.8 Gy/s
	400 V	10 Gy/s	5.0 Gy/s
	500 V	16 Gy/s	7.8 Gy/s
Max. dose per irradiation pulse	300 V	0.69 mGy	0.34 mGy
	400 V	0.91 mGy	0.46 mGy
	500 V	1.14 mGy	0.57 mGy

PTW

IONIZATION CHAMBER TYPE 30004

Radiation Therapy


Diagnostic Radiology

Nuclear Medicine


Radiation Protection

“C/Al-Farmer” Chamber 0.6 cm³

Design




Energy Dependence



Directional Dependence

in air



in PMMA

Reading / Dose [%]

	0°	15°	30°	45°	60°	75°	90°	105°	120°	135°	150°	165°	180°
70 kV	97	97	97	97	97	97	97	97	97	97	97	97	97
100 kV	97	97	97	97	97	97	97	97	97	97	97	97	97
160 kV	97	97	97	97	97	97	97	97	97	97	97	97	97
280 kV	97	97	97	97	97	97	97	97	97	97	97	97	97
Co	97	97	97	97	97	97	97	97	97	97	97	97	97

	0°	15°	30°	45°	60°	75°	90°	105°	120°	135°	150°	165°	180°
70 kV	97	97	97	97	97	97	97	97	97	97	97	97	97
100 kV	97	97	97	97	97	97	97	97	97	97	97	97	97
160 kV	97	97	97	97	97	97	97	97	97	97	97	97	97
280 kV	97	97	97	97	97	97	97	97	97	97	97	97	97
Co	97	97	97	97	97	97	97	97	97	97	97	97	97

in air (without build-up cap) in PMMA (without build-up cap)

PTW

100