



Expect Service

**Radiation Products Design
Inc**

INSTRUCTIONS

RPD INFORMATION

Address	5218 Barthel Industrial Drive Albertville, MN 55301
Website	www.rpdinc.com
Email	sales@rpdinc.com
Phone	763-497-2071 or 800-497-2071
Fax	763-497-2295

RPD PRODUCT INFORMATION

Item Number	Description
970-656-10	10mm Dia. Rounded Nylon Dilator
970-656-15	15mm Dia. Rounded Nylon Dilator
970-656-20	20mm Dia. Rounded Nylon Dilator
970-656-25	25mm Dia. Rounded Nylon Dilator
970-656-30	30mm Dia. Rounded Nylon Dilator
970-656-35	35mm Dia. Rounded Nylon Dilator
970-656-40	40mm Dia. Rounded Nylon Dilator
970-656-45	45mm Dia. Rounded Nylon Dilator
970-656-50	50mm Dia. Rounded Nylon Dilator
946-200	Adjustable Holder/Base for Dilators w/ Stem Scale
946-210	Adjustable Holder and Slessinger Board for Dilators with Stem Scale
978-610	Adjustable Holder/Board for Dilators w/ Stem Scale
946-220	20mm Dia. Rounded Nylon Dilator with Stem Scale
946-225	25mm Dia. Rounded Nylon Dilator with Stem Scale
946-230	30mm Dia. Rounded Nylon Dilator with Stem Scale
946-235	35mm Dia. Rounded Nylon Dilator with Stem Scale

INFO 970-656-20-A

DISCLAIMER

THESE PRODUCTS ARE NOT STERILE AND ARE TO BE USED BY AUTHORIZED PERSONNEL ONLY.

RADIATION PRODUCTS DESIGN INC assumes no liability for consequential damages of any kind for this material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer. Federal law (USA) restricts the sale of this device for use only by (or at the order of) a physician.

USE

Can be used in CT and MR.

CLEANING

Use detergent similar to Klenzyme®, Manu-Klenz®, or NpH Klenz® and water.

PATIENT CLEANING

After use, clean with hot water and soap. Rinse thoroughly.

STERILIZATION

Sterrad	Maximum temperature of unit NOT to exceed 140°F. Type of chemical used for sterilization: Peroxide.
Steris	Maximum temperature of unit NOT to exceed 140°F.
Steam Autoclave	Normal Cycle

End of Document