



Radiation Products Design, Inc.

Item 970-641-2 and 970-641-3 Nylon Ovoid Caps

This Product is not sterile

Nylon Ovoid Cap Fit: The ovoid caps should fit snug on the applicator. If the caps are too loose sterile KY lubricating jelly (or comparable lubricant) may be used inside the cap to provide a suction effect.

Steam Autoclave at Normal Cycle

To be used by authorized personnel only.

Radiation Products Design, Inc. assumes no liability or for consequential damages of any kind for this material or any direct or indirect result and consequences of its use or misuse by the customer.

**5218 Barthel Industrial Drive NE – Albertville, MN 55301
Phone: 763-497-2071 or 800-497-2071 Fax: 763-497-2295**