



Expect Service

Radiation Products Design Inc

INSTRUCTIONS

RPD INFORMATION

Address	5218 Barthel Industrial Drive Albertville, MN 55301
Website	www.rpdinc.com
Email	sales@rpddinc.com
Phone	763-497-2071 or 800-497-2071
Fax	763-497-2295

RPD PRODUCT INFORMATION

Item Number	Description
970-541-50	2.5cm Left/Right Nylon Ovoid Cap
970-541-52	3.0cm Left/Right Nylon Ovoid Cap



DISCLAIMER

THESE PRODUCTS ARE NOT STERILE AND ARE TO BE USED BY AUTHORIZED PERSONNEL ONLY.

RADIATION PRODUCTS DESIGN INC assumes no liability for consequential damages of any kind for this material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer. Federal law (USA) restricts the sale of this device for use only by (or at the order of) a physician.

NYLON OVOID CAP FIT

The ovoid caps should fit snug on the applicator. If the caps are too loose sterile KY lubricating jelly (or comparable lubricant) may be used inside the cap to provide a suction effect.

STERILIZATION

Steam Autoclave at Normal Cycle

End of Document