



Expect Service

Radiation Products Design Inc

INSTRUCTIONS

RPD INFORMATION

Address	5218 Barthel Industrial Drive Albertville, MN 55301
Website	www.rpdinc.com
Email	sales@rpdinc.com
Phone	763-497-2071 or 800-497-2071
Fax	763-497-2295

RPD PRODUCT INFORMATION

Item Number	Description
970-656-10	10mm Diameter Rounded Nylon Dilator
970-656-15	15mm Diameter Rounded Nylon Dilator
970-656-20	20mm Diameter Rounded Nylon Dilator
970-656-25	25mm Diameter Rounded Nylon Dilator
970-656-30	30mm Diameter Rounded Nylon Dilator
970-656-35	35mm Diameter Rounded Nylon Dilator
970-656-40	40mm Diameter Rounded Nylon Dilator
970-656-45	45mm Diameter Rounded Nylon Dilator
970-656-50	50mm Diameter Rounded Nylon Dilator

DISCLAIMER

THESE PRODUCTS ARE NOT STERILE AND ARE TO BE USED BY AUTHORIZED PERSONNEL ONLY.

RADIATION PRODUCTS DESIGN INC assumes no liability for consequential damages of any kind for this material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer. Federal law (USA) restricts the sale of this device for use only by (or at the order of) a physician.

CLEANING

Use detergent similar to Klenzyme®, Manu-Klenz®, or NpH Klenz® and water.

PATIENT CLEANING

After use, clean with hot water and soap. Rinse thoroughly.

STERILIZATION

Sterrad	Maximum temperature of unit NOT to exceed 140°F. Type of chemical used for sterilization: Peroxide.
Steris	Maximum temperature of unit NOT to exceed 140°F.
Steam Autoclave	Normal Cycle

End of Document