

Radiation Products Design, Inc. 5218 Barthel Industrial Drive Albertville, MN 55301 www.rpdinc.com

Phone: 800-497-2071 Fax: 763-497-2295

RPD is an authorized distributor

### **RPD Product Information**

## **Item Number Description**

\* Exradin Model A1 - Miniature Shonka Thimble Chamber 300-205-A1

300-205-A1SL \* Exradin A1SL - Miniature Shonka Thimble Slim Line Chamber

300-205-M1 \* Exradin Model M1 - Miniature Shonka Thimble Chamber

300-205-T1 \* Exradin Model T1 - Miniature Shonka Thimble Chamber

# EXRADIN.


#### IONIZATION

#### **CHAMBER**

**SPECIALISTS** 

# Model 1 - Miniature Shonka Thimble Chamber

- Axially symmetric design
- Small gas volume
- Homogeneous construction
- Complete guarding
- ■Waterproof construction
- ■Gas flow capability
- Integral flexible low-noise cable
- Separable minimum scatter stem
- Matching equilibrium caps
- Handsome wooden case



#### **Features**

Collecting Volume: 0.05 cm<sup>3</sup>

Centroid of Collecting Volume: 3.8 mm from tip of chamber

Collector Diameter: 1.5 mm

Outside Diameter Sensitive Region: 6.0 mm

Wall Thickness: 1.0 mm

Wall, Collector, and Guard Material:

Model A1 – Shonka air-equivalent plastic C552

Model M1 – Magnesium

Model T1 – Shonka tissue-equivalent plastic A150

Maximum Polarizing Potential: Greater than 1000 V

**Inherent Leakage Currents:** Less than 10<sup>-15</sup> Amperes

Cable: 50 Ohms, 29pF/f, 2 m long

Signal Connector: Choice of coaxial and triaxial BNC plugs and jacks and

triaxial TNC plug

**High Voltage Connector:** Integral with triaxial connectors (shell of chamber is common with connector body). Banana, MHV, and SHV plugs are also available.

**Vent Tubes:** Two 21 Ga stainless steel tubes which may be extended to connector via 22 Ga TFE tubing for possible immersion in liquid phantoms.

**Stem:** Black phenolic, 0.5 in OD by 9, 5, or 4 in long. Two-piece stem consisting of a 4 in segment at chamber and a 5 in segment threaded together is also available.