



Expect Service

Radiation Products Design Inc

INSTRUCTIONS

RPD INFORMATION

Address	5218 Barthel Industrial Drive Albertville, MN 55301
Website	www.rpdinc.com
Email	sales@rpdinc.com
Phone	763-497-2071 or 800-497-2071
Fax	763-497-2295

RPD PRODUCT INFORMATION

Item Number	Description
455-300	Cervix 3/16" Stainless Steel Ball Markers with Suture Hole, 50/Pkg
455-301	Cervix 3/16" Delrin Ball Markers with Suture Hole, 50/Pkg

DISCLAIMER

THESE PRODUCTS ARE NOT STERILE AND ARE TO BE USED BY AUTHORIZED PERSONNEL ONLY.

RADIATION PRODUCTS DESIGN INC assumes no liability for consequential damages of any kind for this material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer. Federal law (USA) restricts the sale of this device for use only by (or at the order of) a physician.

PRECAUTION

These products are for single patient use only.

CLEANING

The device should be washed with detergent and water and a small brush should be used to clean out inside the hole prior to sterilization.

STAINLESS STEEL STERILIZATION

Sterilization	Autoclave	Use your manufacturer's recommendations or procedures for stainless steel surgical instruments.
	Sterrad NX100	Use your manufacturer's recommendations or procedures for stainless steel surgical instruments.

DELTRIN STERILIZATION

Sterilization	Autoclave	Use your manufacturer's recommendations or procedures for delrin surgical instruments.
	Sterrad NX100	Use your manufacturer's recommendations or procedures for delrin surgical instruments.