



Expect Service

Radiation Products Design Inc

INSTRUCTIONS

RPD INFORMATION

Address	5218 Barthel Industrial Drive Albertville, MN 55301
Website	www.rpdinc.com
Email	sales@rpdinc.com
Phone	763-497-2071 or 800-497-2071
Fax	763-497-2295

RPD PRODUCT INFORMATION

Item Number	Description
495-604	Wedge Tray Scale f/ TBI Isodose Curve, Varian III

INTRODUCTION

The Radiation Physicist must find the projected isodose curves for their TBI Stand location.

The mapping of the isodose curves must be marked on the wall or stand covered with paper so values and lines can be seen. Use chamber or film. Record gantry angle used.

The Total Body Scale Tray is inserted in the treatment machine and rotated to the same gantry angle setting used for producing the projected isodose on the wall or stand. With a black marker, mark the demagnification on the acrylic TBI tray.

End of Document