



Diamond Clamp

Item #	Clamp Width	Hole A Diameter	Hole B Diameter	Hole C Diameter
16	.49" (12.5mm)	3.5-5mm	6-8mm	9-12mm
17	.49" (12.5mm)	9-12mm	.6-8mm	3.5-5mm




One Hole

Item #	Clamp Width	Hole A Diameter	Distance to Tip
05	.49" (12.5mm)	8mm	.182" (4.6mm)
12	.49" (12.5mm)	5mm	.182" (4.6mm)




Two Holes

Item #	Clamp Width	Hole A Diameter	Distance to Tip	Hole B Diameter	Distance to Tip
01	.49" (12.5mm)	6mm	.182" (4.6mm)	8mm	.579" (14.7mm)
02	.49" (12.5mm)	8mm	.182" (4.6mm)	6mm	.579" (14.7mm)
03	.49" (12.5mm)	6mm	.182" (4.6mm)	12mm	.657" (16.7mm)
04	.49" (12.5mm)	12mm	.306" (7.8mm)	6mm	.781" (19.8mm)
06	.25" (6.4mm)	.25" (6.4mm)	.226" (5.8mm)	.312" (7.9mm)	.632" (16.1mm)
07	.25" (6.4mm)	.312" (7.9mm)	.226" (5.8mm)	.25" (6.4mm)	.632" (16.1mm)
10	.49" (12.5mm)	6.3mm	.149" (3.8mm)	5.8mm	.512" (13.0mm)
14	.49" (12.5mm)	.25" (6.4mm)	.182" (4.6mm)	.375" (9.5mm)	.549" (14.0mm)



Three Holes

Item #	Clamp Width	Hole A Diameter	Distance to Tip	Hole B Diameter	Distance to Tip	Hole C Diameter	Distance to Tip
08	.49" (12.5mm)	4mm	.118" (3mm)	6mm	.377" (9.6mm)	9mm	.735" (18.7mm)
15	.49" (12.5mm)	6mm	-	8mm	-	10mm	-



One Hole & One Slot

Item #	Clamp Width	Hole A Diameter	Distance to Tip	Slot Length and Diameter	Distance to Tip
09	.49" (12.5mm)	6.5mm	.165" (4.20mm)	18mm x 8mm	1.11" (28.2mm)
11	.49" (12.5mm)	6.3mm	.182" (4.63mm)	6mm x 5mm	.549" (14.0mm)


How to Order

Order the arm system with the item number extension added to the end of the item number for the clamp that matches your intended use

Example: Choose item number 978-60501 for MR/CT Titanium Flex Arm System with Wide 6mm/8mm Clamp Holes


Note: Custom clamps can be made to order if your desired specifications are not listed above.





Special Fixation Element Clamps

Item #	Clamp Width	Fixation Element
13	.49" (12.5mm)	Elekta Fixation Element #126024
18	.49" (12.5mm)	Elekta Venezia™ Advanced Gynecological Applicator



Special Fixation Elements

These clamps are specifically designed to work with other established 3rd party GYN systems.

How to Order

Order the arm system with the item number extension added to the end of the item number for the clamp that matches your intended use

Example: Choose item number 978-60501 for MR/CT Titanium Flex Arm System with Wide 6mm/8mm Clamp Holes

Note: Custom clamps can be made to order if your desired specifications are not listed above.



PHONE (763) 497-2071



SALES@RPDINC.COM



FAX (763) 497-2295



WWW.RPDINC.COM



5218 BARTHEL INDUSTRIAL DRIVE
ALBERTVILLE, MN 55313