



**Expect Service**

**Radiation Products Design Inc**

## **INSTRUCTIONS**

### **RPD INFORMATION**

<b>Address</b>	5218 Barthel Industrial Drive Albertville, MN 55301
<b>Website</b>	<a href="http://www.rpdinc.com">www.rpdinc.com</a>
<b>Email</b>	<a href="mailto:sales@rpdinc.com">sales@rpdinc.com</a>
<b>Phone</b>	763-497-2071 or 800-497-2071
<b>Fax</b>	763-497-2295

### **RPD PRODUCT INFORMATION**

<b>Item Number</b>	<b>Description</b>
067-650	DOT Source Transport Pig, RENTAL ONLY

### **DISCLAIMER**

**THESE PRODUCTS ARE NOT STERILE AND ARE TO BE USED BY  
AUTHORIZED PERSONNEL ONLY.**

RADIATION PRODUCTS DESIGN INC assumes no liability for consequential damages of any kind for this material when used interchangeably with products of other manufacturers/suppliers or for any direct or indirect results and consequences of its use or misuse by the customer. Federal law (USA) restricts the sale of this device for use only by (or at the order of) a physician.

## FULL ASSEMBLY



## INSTRUCTIONS FOR USE



Step 1: Remove Cover



Step 2: Turn Container Upside Down



Step 3: Remove Container



Step 4: Remove Top of Lead Pig