

## ALRA® RADIATION THERAPY LOTION



- Provides quick, reliable relief from acute radiation dermatitis
- Non-irritating, non-allergenic
- Does not interfere with ongoing treatment
- Excellent patient acceptance

Alra® Radiation Therapy Lotion is specifically developed for the treatment of skin exposed to radiation therapy. It contains a highly concentrated emulsion of aloe vera gel, lanolin, vitamin E, soluble collagen, and allantoin. It does not contain water, glycerin, mineral oil, alcohol or added fragrance.

In patient studies conducted at the Cedars-Sinai Comprehensive Cancer Center in Los Angeles, none of the patients treated with Alra® Radiation Therapy Lotion showed irritation or allergic reactions.

This soothing lotion aids in the replacement of damaged tissue, helping the natural healing process to take place.

Item #	Description	Quantity
125-000	Alra® Radiation Therapy Lotion, 4 oz. Bottle	1
125-001	Alra® Radiation Therapy Lotion, 4 oz. Bottle	12
125-002	Alra® Radiation Therapy Lotion, 1 oz. Bottle	1
125-003	Alra® Radiation Therapy Lotion, 1 oz. Bottle	24

## ALRA® NON-METALLIC DEODORANT



- Free of metallic ingredients that can interfere with radiation therapy treatments

Alra® Non-Metallic Deodorant is specifically formulated to eliminate perspiration odor and wetness without the use of metallic ingredients that may interfere with treatments. This gentle, non-irritating deodorant contains oil of cypress, a clean-smelling natural fragrance, and lichen extract which is a natural bactericide and fungicide.

Item #	Description	Quantity
125-005	Alra® Non-Metallic Deodorant	1
125-006	Alra® Non-Metallic Deodorant	12

## MEDI-SOOTHE™ LOTION with 80% “Diamond-Pure” Aloe Vera



- 80% “Diamond-Pure” Aloe Vera
- Vitamins E and B-5
- Alcohol Free
- Non greasy
- Allantoin promotes healing, speeds cell regeneration
- Non-prescription formula does not contain waxes or metals

Medi-Soothe™ lotion provides soothing relief when used to help heal various forms of skin breakdown; Radiation Induced Dermatitis, thermal burns, 1st & 2nd degree burns, sunburn and rashes.

### Benefits

- Effective barrier to heat build-up
- Replenishes lost moisture
- Encourages production of collagen fibrillation
- Speeds cell regeneration
- Overall promotes healing

Item #	Description	Size
126-002	Medi-Soothe™ Lotion	2 oz Bottle
126-006	Medi-Soothe™ Lotion	6 oz Flip Top Bottle

## ALRA® MILD CONDITIONING SHAMPOO



- Minimizes hair damage while cleansing and conditioning hair and scalp

Oncology therapy can cause normal hair to become extremely dry and brittle, and the scalp to be very sensitive. While no shampoo can prevent hair loss caused by some chemotherapy drugs, Alra® Mild Conditioning Shampoo is an extra-gentle formula that is acid-balanced and specially-formulated to minimize hair damage.

Ingredients include: aloe vera gel, DL panthenol (pro-vitamin B-5) and aromatic oil of rosemary, to stimulate scalp circulation.

Item #	Description	Quantity
125-010	Alra® Mild Total Body Shampoo, 8 oz Bottle	1
125-011	Alra® Mild Total Body Shampoo, 8 oz Bottle	12

# PATIENT & PERSONAL PRODUCTS

## RADIA-GUARD® METAL-FREE RADIATION TREATMENT AIDS



**Item 124-000 Radia-Guard® Discreet Deodorant** is aluminum and metallic free. Discreet spray deodorant is made with all-natural ingredients while still giving you the odor-fighting effectiveness you need. Rest assured- it is completely free from aluminum and other metallic compounds. Can be applied to underarms as often as needed to control odor. Avoid spraying on open wounds. *Ingredients: SDA 40 alcohol, deionized water, polysorbate 20, citrus lemon oil, lavender oil, sage oil.*

your skin. Jojoba and other moisturizers work to help keep your skin healthy and deeply moisturized, which can help with skin elasticity after your treatment cycle is completed. Radia-Guard® Treatment Lotion should be used from the beginning of your treatment cycle for best results and can be applied immediately after every treatment. Use more frequently if skin becomes irritated due to your treatments. This lotion absorbs in minutes and does not leave a greasy or sticky residue. Avoid applying the lotion on any temporary skin markings. If severe skin irritation or open sores develop as a result of therapy discontinue use and consult your doctor. *Ingredients: DI water, aloe vera gel, emulsifying wax, jojoba oil, cyclomethicone, lidocaine HCl, methyl paraben, ethyl paraben, butyl paraben, almond oil, cetal, carbomer, triethanolamine, green tea extract (leaf).*

**Item 124-005 Radia-Guard® ReStructure Cream** is a unique blend of the richest moisturizers and other choice ingredients that provide healing and revitalizing benefits to skin damaged by radiation. ReStructure Cream will help regain elasticity, improve skin tone and smoothness, and protect skin from damaging free radicals. DMAE tightens pores and fights wrinkles. *Ingredients: Purified water, glyceryl mono stearate, jojoba oil, shea butter, methyl gluceth 20, methyl glucose sesquisteerate, Peg-20 glucose sesquisterate, cetyl acetate, acetylated lanolin alcohol, polysorbate 80, royal jelly, dimethylaminoethanol (DMAE), methyl paraben, fragrance, propyl paraben, DMDM hydantoin, vitamin E.*

**Item 124-015 Radia-Guard® Heavenly Lip Balm** is made with the best natural ingredients and is perfect for relieving the dry chapped lips of radiation therapy patients. *Ingredients: sweet almond oil, coconut oil, beeswax, castor oil, unrefined cocoa butter, mango butter, vitamin E, comfrey oil, and flavoring.*

**Item 124-010 Radia-Guard® Radiation Treatment Lotion** is completely metal-free and contains Aloe and Licocaine to soothe

Item #	Radia-Guard®	Quantity
124-000	Discreet Deodorant	12
124-005	ReStructure Cream	12
124-010	Radiation Treatment Lotion	12
124-015	Heavenly Lip Balm	12

## SPAND-GEL™

### Sterile Hydrogel Sheets with Diamond Aloe Vera



Spand-Gel Hydrogel sheets are comprised of cross-linked water swollen hydrogel polymers which contain deionized water, polyethylene oxide and aloe vera.

#### Hydrogel sheets are an effective, closed-dressing procedure

- It's highly conformable and permeable
- Absorbs varying amounts of drainage, depending on composite
- Spandage™ effectively holds the face mask in place

#### Spand-Gel™ Sterile Primary Wound Dressing provides:

- A moist, cool wound-healing environment
- Promotes tissue granulation formation
- Relieves pain, burning, itching and soreness
- Absorbs wound exudate 3 times its weight
- Provides a barrier from exogenous bacteria and fluids
- Does not adhere to the tissue



- Full and half face masks sized for entire treatment area
- Neck wrap is contoured and reversible
- Clear film on both sides for left or right placement

**Contraindication:** Infected wounds, unless dressing is used in conjunction with appropriate topical antibiotic.

Item #	Spand-Gel™	Quantity
126-010	Sterile 3" x 8" Sheet	30
126-011	Sterile 4" x 4" Sheet	30
126-012	Sterile 5" x 9" Sheet	30
126-013	Sterile 6" x 8" Sheet	30

Item #	Spand-Gel™	Quantity
126-014	Sterile Full Face Mask	6
126-015	Sterile Half Face Mask	6
126-016	Sterile Breast Wrap	6
126-017	Sterile Neck Wrap	6

# PATIENT & PERSONNEL PRODUCTS

## DISPOSABLE VINYL EXAM GLOVES



Tough 5-mil Vinyl Exam Gloves fit either hand. 11" length. Comes lightly powdered or powder free. Excellent general-purpose disposable gloves. Sold in 100/box.

## DISPOSABLE NITRILE GLOVES



These patented gloves provide three times the puncture resistance of natural rubber and vinyl disposable gloves of comparable thickness. Gloves contain no natural rubber latex proteins. Disposable, 100% nitrile gloves are 4 mil thick with low-powder, 9 1/2" long, fits either hand. Sold in 100/box.

Item #	Gloves, Disposable Vinyl	Size
280-010	Powdered	Medium
280-011	Powdered	Large
280-012	Powdered	X-Large
280-015	Powder-Free	Medium
280-016	Powder-Free	Large
280-017	Powder-Free	X-Large

Item #	Gloves, Disposable Nitrile	Size
280-001	Powdered	Medium
280-002	Powdered	Large
280-003	Powdered	X-Large
280-004	Powder-Free	Medium
280-005	Powder-Free	Large
280-006	Powder-Free	X-Large

## LATEX SURGICAL GLOVES



Available powder-free or with an absorbable dusting power that meets USP requirements. Micro-textured. Feature beaded cuffs to prevent rolldown. Right and left hand specific sizes for better fit. Test levels exceed ASTM D 357, ASTM D 5151 (AQL 1.0), ASTM D 6184, and ASTM D 5712. Single use. Sterilized by gamma radiation. Length: 12". 8.8 mil thick. Sold in 50 pairs/pack.

Item #	Gloves, Latex Surgical	Size
280-051	Powdered	6
280-052	Powdered	6 1/2
280-053	Powdered	7
280-054	Powdered	7 1/2
280-055	Powdered	8
280-056	Powdered	8 1/2
280-057	Powdered	9
280-061	Powder-Free	6
280-063	Powder-Free	7
280-064	Powder-Free	7 1/2
280-065	Powder-Free	8
280-066	Powder-Free	8 1/2
280-067	Powder-Free	9

## BIOGEL LATEX SURGICAL GLOVES Powder-Free



With patented Biogel coating. 100% physically inspected. Easy to don, even over damp hands. Curved finger design helps prevent hand fatigue. Beaded cuffs. Hand specific. Single-use. Natural color. Sterile. 10 1/2" length. 10 mil thick. Specially processed to ensure ≤50µg of water-extractable protein per gram of latex. Sold in 50 pairs/pack except for size 8 1/2" which is only available in a case of 200 pairs.

Item #	Gloves, Biogel Latex Surgical	Size
280-071	Powder-Free	6
280-072	Powder-Free	6 1/2
280-073	Powder-Free	7
280-074	Powder-Free	7 1/2
280-075	Powder-Free	8
280-076	Powder-Free	8 1/2

## SECA 216 WALL MOUNTED STADIOMETER

The seca 216 Wall-Mounted Stadiometer is equally suitable for measuring both children and adults. The plastic frame is simply mounted on the wall at the desired height. The stadiometer has a removable scale with a measuring range of 0 to 230 centimetres. The scale can be inserted variably into a 150 centimeter long plastic profile. The user decides which measuring range appears in the profile.

The results are clearly displayed in a window on the side of the rod and held in place by the locking screw on the head piece. Suitable for international use, the scale is printed in centimeters on one side of the tape and centimeters and inches on the other.

### Specifications

**Measuring Range:** Max. 60" / 150 cm variably selectable from 4 - 87" (10 - 220 cm)


**Graduations:** 1/16" / 1mm

**Dimensions:** 4.7" W x 60" H x 8.5" D (12 x 150 x 21.7 cm)

**Weight:** 4.4 lb (2 kg)

CE 0123

Item #	Description
100-216	SECA 216 Wall Mounted Stadiometer
100-217	Replacement Head Piece for Item 100-216



Y

## SECA 222 TELESCOPIC MEASURING ROD

### Wall-mounted with extra wide slide

Specially developed for wall mounting, these aluminum measuring rods are aligned to work exactly to the millimeter from head to toe. The broad measuring slide and heel positioner are professional prerequisites. The telescopic measuring rod ensures that the result can be easily read at eye level even for the tallest patient. When not in use, the measuring slide can be folded down for safety.

### Specifications

**Functions:** Wall fastening essential/possible

**Measuring Range:** 2" - 90.5" (6 - 230 cm)

**Graduation Length:** 1/8" (1 mm)

#### Dimensions

**Unfolded:** 12.1" W x 91.3" H x 11.3" D (30.8 x 232 x 28.7 cm)

**Folded:** 12.1" W x 52" H x 11.3" D (30.8 x 132 x 28.7 cm)

**Weight:** 2.4 lb (1.1 kg)

CE 0123

Item #	Description
100-222	SECA 222 Telescopic Measuring Rod, Inches
100-223	SECA 222 Telescopic Measuring Rod, cm



## SECA 700 MECHANICAL COLUMN SCALE



- Extra high capacity of 500 lb (220 kg)
- Large, low platform
- Simple to operate and easy to read due to eye-level display
- Sturdy construction with cast-iron base
- Integrated Transport castors for mobile use

With its capacity of 500 lb (220 kg), the SECA 700 has been specially designed for weighing overweight and obese patients. To ensure durability, special attention was paid to the high-quality construction. The generously sized, low level platform makes the SECA 700 easy to mount. The weight scale, has been affixed at eye level, is printed on both sides and has large numerals so that the result is easier to read. A special coating reliably prevents the scale from fading even after years of use. The base of the SECA 700 is made of cast iron and has a special coating on the lever apparatus which protects the inner parts against corrosion. The scale can be moved easily and quickly on the integrated transport castors to wherever it is needed. When fitted with an optional SECA 220 Telescopic Measuring Rod (Item 100-265 & 100-266), the SECA 700 makes time-saving weighing and measuring possible in one step.

### Specifications

#### Items 100-250 and 100-251 SECA 700 Mechanical Column Scale

**Functions:** Reset-to-zero function, optional SECA 220 Telescopic Measuring Rod

**Capacity:** 500 lb (220 kg)

**Graduation:** 1/8 lb (50 g)

**Dimensions:** 20.5" W x 61.3" H x 20.5" D (52.1 x 155.7 x 52.1 cm)

**Platform:** 13.2" W x 3.1" H x 13.6" D (33.5 x 7.9 x 34.5 cm)

**Weight:** 35.7 lb (16.2 kg)

#### Items 100-265 and 100-266 Optional

#### SECA 220 Telescopic Measuring Rod

**Measuring Range:** 24 - 78" (60 - 200 cm)

**Graduation Length:** 1/8" (1 mm)

**Dimensions:** 2.2" W x 58.3" H x 11.3" D

(5.6 x 148.1 x 28.7 cm)

**Weight:** 1.5 lb (0.7 kg)

CE0123

Item #	Description
100-250	SECA 700 Mechanical Column Scale, lbs
100-251	SECA 700 Mechanical Column Scale, lbs/kg
100-265	SECA 220 Telescopic Measuring Rod, cm
100-266	SECA 220 Telescopic Measuring Rod, inches

## SECA 769 DIGITAL COLUMN SCALE



- Site-independent as battery operated
- Easy to transport
- Cost-effective and low-maintenance due to low power consumption
- Intelligent functions for versatile use

The SECA 769 Digital Column Scale with BMI Function is equipped with several intelligent functions. The BMI function permits a reliable evaluation of the nutritional condition of the patient. The auto-HOLD function ensures that the result is still shown on the display after the patient has stepped down from the scale. Additional items which are not to be included in the weight result are simply deducted by means of the TARE function. The SECA 769 is extremely energy-saving and thus very cost-effective. The SECA 769 is battery-operated, no wall outlet connection is necessary. The automatic switch-off function is another reason for the low energy consumption. The column scale is supplied with a measuring rod so that both weight and height can be determined in just one time-saving step. With its smooth running transport castors, the column scale is also easy to move around.

### Specifications

#### Item 100-261 SECA 769 Digital Column Scale with BMI Function

**Capacity:** 450 lb (200 kg)

**Graduation:** 0.2 lb (100 g)

**Dimensions:** 11.6" W x 32.7" H x 16.4" D (29.4 X 83.1 X 41.7 cm)

**Platform:** 10.7" W x 3" H x 11" D (27.2 x 7.5 x 28 cm)

**Measuring Range:** 24 - 78" (60 - 200 cm)

**Graduation Length:** Measuring rod 1/8" (1 mm)

**Weight:** 13.7 lb (6.2 kg)

**Power Supply:** Batteries, power adapter optional

**Functions/Accessories:** TARE, auto-HOLD, BMI, lbs/kg switch-over, automatic switch-off, measuring rod

**Optional:** Switch-mode power adapter SECA 400

#### Item 100-263 Optional A/C UL Adapter

**Dimensions:** 1.6" W x 1.4" H x 2.9" D

(4.2 x 3.5 x 7.3 cm)

**Cable length:** 72.8" (1.85 m)

**Weight:** 0.34 lb (0.15 kg)

CE0123

Item #	Description
100-261	SECA 769 Digital Column Scale with BMI Function
100-263	A/C UL Adapter for Item 100-261

# PATIENT & PERSONAL PRODUCTS

## Blood Pressure Monitor



The professional stethoscope can be used for other diagnostic procedures. The self-taking method is easily adapted with the D-Ring cuff. Arm size: 9" to 13".

### Item 288-200 Blood Pressure Monitor Includes

- Gauge
- D-ring cotton cuff with chest piece
- Bladder
- Bulb and air release valve
- Binaural and eartips
- Y-tubing
- Vinyl storage case
- User manual
- Record chart
- 1 year warranty (limited, see manual for details)

Item #	Description
280-200	Blood Pressure Monitor

## FLEXBAK®



Flexbak® helps reduce back stress and strain with a design so unique its patented. Its three-part support is the key to a system that provides maximum back support with ease of movement. The 9" medical-style truss uses two internal support panels above and below the waist surrounded by an elastic third panel. These combine to produce maximum bracing without pinching or riding up.

To determine size, measure around abdomen at the navel.

Item #	Description	Size
280-11501	Flexbak® - Small	26" - 36"

## FLEXRIST®



FlexRist®, the comfortable all around support made of Neoprene is the best defense against Carpal Tunnel Syndrome. One design can be adapted to any wrist, right or left. It washes easily, is comfortable, extremely lightweight and versatile.

Item #	Description
280-11101	FlexRist® Regular

## MAXRIST®



MaxRist® offers a restrictive metal spoon stay that fits in the palm to eliminate the symptoms of CTS and other repetitive strain injuries. The lightweight flexible material does not hinder finger movement.

To determine size, measure around the wrist for correct size.

Item #	Description	Size
280-10902	MaxRist® Right	Medium (6 1/2" to 7 1/2")
280-11002	MaxRist® Left	Medium (6 1/2" to 7 1/2")

## RIST RITE®



The Rist Rite® is a neoprene support that is reinforced with Cordura nylon across the palm for long wear. It has one strap that stops short of the hand.

To determine size, measure around the wrist.

Item #	Description	Size
280-11202	Rist Rite® Right	Medium (6 1/2" to 7 1/2")
280-11203	Rist Rite® Left	Medium (6 1/2" to 7 1/2")